



## COVID-19 Antigen Home Test

A rapid test for the detection of COVID-19 antigens in nasal specimens.







## COVID-19 Antigen Home Test


A simple way to test of COVID-19

# Get Back to Sharing Special Moments *with Family and Friends*

The Flowflex COVID-19 Antigen Home Test is all you need to determine your family's COVID-19 status.


### Find peace of mind with Flowflex!

- Easy-to-use nasal swab test
- Can be used to test children as young as 2 years old
- For use with and without COVID-19 symptoms
- Results in 15 minutes
- No need to send off to a lab to obtain results
- Compact packaging for "On-The-Go" testing







## Test Procedure Overview



**STEP 1**  
**Swab**



**STEP 2**  
**Swirl**



**STEP 3**  
**Drop**




**STEP 4**  
**See Results**






# Flowflex<sup>®</sup>

**1PK TEST**  
**NDC NUMBER:**  
82607-0660-26 (1pk)



Fast



Easy to Use



Trusted



KYC LLC, 1000 WILKINSON TRACE, BOWLING GREEN, KY 42103 | 270-407-5400 | WWW.KYC-LLC.COM



**KYC LLC, 1000 WILKINSON TRACE, BOWLING GREEN, KY 42103 | 270-407-5400 | WWW.KYC-LLC.COM**